




ResusRight
FOR SAFER BIRTH

Skills Assessment Mode



Short press on the ⓘ button opens the menu of the JUNO. Press the right arrow button to navigate to **Skills Assessment**, short press on the ⓘ button to select.




JUNO now shows the manikin size selected that will be used to calculate effectiveness of ventilation. Change if necessary.




Follow the instructions on the screen. The assessment will start when ventilating. During the assessment you will not receive feedback from the JUNO. When finished press the ⓘ button.


How to interpret your results




The JUNO will show how many ventilations were in the target range for each parameter for the manikin size.



When continuing to the next screen the JUNO will show you more detail on the Mask Seal.



The next screen will show you information on the Tidal Volumes delivered.



On the final screen details on the Respiratory Rate are shown.

When pressing finish, the JUNO will go back to the start of Skills Assessment mode.



Scan QR code or go to www.resusright.com/support for troubleshooting.

Please contact us at support@resusright.com for further support.

Device Information



The JUNO Training Monitor is designed to improve the effectiveness of skills training on neonatal manikins. Current neonatal manikins provide no objective feedback around mask technique or volume delivery, making it hard to effectively teach skills and ensure every clinician is ready for a newborn resuscitation. The JUNO clips between the face mask and the resuscitative device to provide feedback on correct resuscitation technique, assisting educators to improve the quality of training on neonatal manikins.

WARNING: DEVICE IS NOT TO BE USED ON HUMANS, FOR TRAINING PURPOSES ONLY

Parameters	Accuracy
Tidal Volume	$\pm 0.5\text{mL}$ or $\pm 8\%$
Respiratory Rate	± 1 bpm

ResusRight Pty Ltd
Sydney Knowledge Hub
Level 2, Mereweather Building H04,
The University of Sydney NSW 2006
www.resusright.com

Setting up the JUNO



Ensure the JUNO is charged prior to use. Unplug the charger before using the JUNO.



Select the correct size mask for your training manikin and connect the mask to the JUNO sensor.


When using a triangular mask, ensure the tip of the mask is pointing towards the device.



Attach the T-piece or self-inflating bag to the other side of the sensor. Adjust so the resuscitation device can be comfortably operated.

Training Mode



Press and hold the  button to power on the device.
The device will show a blue ResusRight screen when powering up, you can now release the power button.



The JUNO is ready for use when the **Training Screen** is shown.



Apply the mask over nose and mouth and start ventilating as normal. The JUNO will provide breath feedback as shown on next page.

Training Mode Information

LOW LEAK

2 - 3 kg

VTE **5.2** mL

RATE **60** bpm

N/A

NO LEAK

LOW LEAK

HIGH LEAK

No ventilation detected

Less than 30% leak

Between 30 - 60% leak

More than 60% leak

Indicates the size of the infant that the delivered tidal volumes are appropriate for.

VTE **5.2** mL

Vte in normal range

RATE **60** bpm

Ventilatory rate below 80bpm

VTE **2.2** mL

Turns red if Vte is too high or low

RATE **200** bpm

Turns red if ventilatory rate is above 80bpm